



Ref. 1469





GENERAL CHARACTERISTICS

Is modern, desktop laboratory centrifuge. Its structure guarantees easy service, safe work and the wide range of applications in medical, biochemical and others laboratories, where constituents of investigated material are to be separated under the effect of centrifugal force.

- Microprocessor control of spinning parameters provides high repeatability of results,
- LCD display of actual parameters: speed, time and relative centrifugal force. Optimization of acceleration and deceleration force.
- “Short” spinning function.
- Unbalance protection system.
- It is not possible to open the lid, during spinning.
- It is not possible to start, when the lid is open.
- Full metal case and stainless steel chamber.

Product according to national and international safety regulations EN-61010-1 and EN-61010-2-020 for laboratory equipments (CE).

Ref	Dimensions W x H x D	Weight	Consumption W
1469	410x370x510mm.	35Kg	225
Capacity	Max. R.p.m	Max. RCF	
1000 (4x250) ml.	4300	3013 xg.	
Maintenance-free	Timer	Voltage	Brushless
Yes	1-99 min ∞	130/250v. 50/60Hz.	Yes
Display	Rotor identification	Refrigerant	Temperature range
LCD	Yes	NO	---

Rotor		Accessories information			
 <p>Ref. 1469-1253C</p>	24x15ml	Accessories	Ø x Long	Tubes	Ref.
	Angle 30 Max. 4300 R.p.m. Max. 2233 xg	Round carriers 5ml.	Ø13x75mm.	24	4088
		Round carriers 7ml.	Ø13x100mm	24	4082
Rotor		Accessories information			
 <p>Ref. 1469-1254</p>	10x50ml Fal.	Accessories	Ø x Long	Tubes	Ref.
	Angle 30 Max. 4300 R.p.m. Max. 2749 xg	Buckets not included PC cup.	Ø30x100mm	10	3275C
		Buckets not included	Ø30x100mm	10	3276
	Round carriers 30 / 25ml.	Ø 25x100mm.	10	4248	
Rotor		Accessories information			
 <p>Ref. 1469-2436a</p>	4x200ml.	Accessories	Ø x Long	Tubes	Ref.
	Swing out Max. 4300 R.p.m. Max. 3013 xg	Buckets for 2x50ml. Fal.	Ø30x120mm	8	3042
		Buckets for 50ml. Fal.	Ø30x120mm	4	3045
		Hangers with buckets for 15ml. Fal	Ø16x120mm	16	3044C
		Hangers with buckets for 15ml. Fal	Ø16x120mm	8	3129C
		Buckets 200ml.	Ø56.5x94mm	4	3437
		Buckets W/C 200ml.	Ø56.5x94mm	4	3438C
		Round carriers 100ml.	Ø45.5x100mm	4	4104
		Round carriers 5 / 7ml.	Ø13x75/100mm	28	4106
		Round carriers 5 / 7ml.	Ø13x75/100mm	32	4107
		Round carriers 10ml.	Ø16.7x75mm	28	4108
		Round carriers 5ml.	Ø13.5x85mm	28	4109
		Round carriers 15 / 10ml.	Ø16.7x100mm	28	4110
		Round carriers 15 ml.	Ø16.7x110mm	20	4111
		Round carriers 30 / 25ml.	Ø25x100mm	12	4112
		Round carriers 50ml. Fal.	Ø30x120mm	4	4113
		Round carriers 7ml. / without cap	Ø12.1x100mm	48	4441
		Round carriers 50ml.	Ø35x100mm	4	4444
		Round carriers 5ml. / without cap	Ø12.1x75mm	48	4446
		Round carriers 12 / 15ml.	Ø17/19.4x107mm	16	4449
Round carriers 2.2 / 1.5ml.	Ø11x38.5mm	36	4450		
Rotor		Accessories information			
 <p>Ref. 1469-2451C</p>	2x3 microtiter	Accessories	Ø x Long	Tubes	Ref.
	Swing out Max. 3000 R.p.m. Max. 1467 xg				