



Ref. 1463 R

GENERAL CHARACTERISTICS.

GOLGI-R is modern, desktop, refrigerated laboratory microcentrifuge. Its structure guarantees easy service, safe work and the wide range of applications in medical, biochemical and other laboratories, where constituents of investigated material are to be separated under the effect of centrifugal force.

- Microprocessor control of spinning parameters provides high repeatability of results,
- LCD display of actual parameters,
- Free programming of speed or relative centrifugal force (RFC), time, temperature, characteristics of acceleration and deceleration,
- Spinning in "short" mode.
- Inspection glass.
- Errors indication system.
- Low interference level (Power Factor Module).
- Ecological and efficient cooling system.
- Pre-freezing mode provides stabilization of temperature before and after spinning, Compressor overloading protection system.

Product according to national and international safety regulations EN-61010-1 and EN-61010-2-020.

Ref	Dimensions W x H x D	Weight	Consumption W
1463 R	280x298x564mm.	35Kg	550
Capacity	Max. R.p.m	Max. RCF	
79.2 ml.	300-14000	18188 xg.	
Maintenance-free	Timer	Voltage	Brushless
Yes	1-99 min ∞	115/230v. 50/60Hz.	Yes
Display	Rotor identification	Refrigerant	Temperature range
LCD	Yes	Yes	-20 + 40

Rotor		Accessories information			
 <p>Ref. 1463 R-1731/1732</p>	24x1,5/2,2ml.	Accessories	Ø x Long	Tubes	Ref.
	36x1,5/2,2ml			24 / 36	
	Angle 45 Max. 14000 R.p.m. Max. 18188 xg Hermetically sealed lid	Round carriers 0,2 ml.	Ø6,2mm	24 / 36	4133
		Round carriers 0,4 ml.	Ø5,7mm	24 / 36	4126
		Round carriers 0,5 ml.	Ø8mm	24 / 36	4084
Rotor		Accessories information			
 <p>Ref. 1463 R-1738/1735</p>	18x1,5/2,2ml	Accessories	Ø x Long	Tubes	Ref.
	24x1,5/2,2ml			18 / 24	
	Filter				
	Angle 45 Max. 14000 R.p.m. Max. 18188 xg Hermetically sealed lid				
Rotor		Accessories information			
 <p>Ref. 1463 R-1733</p>	36x0,5ml.	Accessories	Ø x Long	Tubes	Ref.
	Angle 45 Max. 14000 R.p.m. Max. 17968 xg Hermetically sealed lid			36	
		Round carriers 0,2 ml.	Ø6,2mm	36	4133
Rotor		Accessories information			
 <p>Ref. 1463 R-1734/2730C</p>	4x(8x0,2ml.) PCR	Accessories	Ø x Long	Tubes	Ref.
	24x1,5/2,2ml			32/24	
	Angle 45 Swing out Max. 14000 R.p.m. Max. 15339/15996 xg Hermetically sealed lid				
Rotor		Accessories information			
 <p>Ref. 1463-2480</p>	24 capilares.	Accessories	Ø x Long	Tubes	Ref.
	Ø1,4x75mm		Ø1,4x75mm	24	
	Horizontal Max. 13000 R.p.m. Max. 16900 xg Hermetically sealed lid	Flat hematocrite reader included			6135
		Round hematocrite reader			6150