

Rotor		Accessories information			
	24x1,5/2,2ml.	Accessories	Ø x Long	Tubes	Ref.
	Angle 45 Max. 14000 R.p.m. Max. 18188 xg Hermetically sealed lid	Round carriers 0,2 ml.	Ø6,2mm	24	4133
		Round carriers 0,4 ml.	Ø5,7mm	24	4126
		Round carriers 0,5 ml.	Ø8mm	24	4084
Ref. 1463-1731					
Rotor		Accessories information			
	18x1,5/2,2ml	Accessories	Ø x Long	Tubes	Ref.
	Filter			18	
	Angle 45 Max. 14000 R.p.m. Max. 18188 xg Hermetically sealed lid				
Ref. 1463-1738					
Rotor		Accessories information			
	36x0,5ml.	Accessories	Ø x Long	Tubes	Ref.
	Angle 45 Max. 14000 R.p.m. Max. 17968 xg Hermetically sealed lid			36	
		Round carriers 0,2 ml.	Ø6,2mm	36	4133
Ref. 1463-1733					
Rotor		Accessories information			
	4 x (8x0,2ml.)	Accessories	Ø x Long	Tubes	Ref.
	PCR			32	
	Angle 45 Max. 14000 R.p.m. Max. 15339 xg Hermetically sealed lid				
Ref. 1463-1734					
Rotor		Accessories information			
	24 capilares.	Accessories	Ø x Long	Tubes	Ref.
	Ø1,4x75mm		Ø1,4x75mm	24	
	Horizontal Max. 13000 R.p.m. Max. 16900 xg	Flat hematocrite reader included			6135
		Round hematocrite reader			6150
Ref. 1463-2480					