



Ref. 1460






## GENERAL CHARACTERISTICS

Is modern, médium capacity centrifuge. It's intended for laboratories where there aren't a big number of samples. It has a wide range of exchangeable rotors ( Capacity on angle rotor 10x50ml Falcon and 24x15ml; swing out rotor from 4x100ml to 48x5/10ml) and wide range of other accessories.




- Microprocessor control of spinning parameters provides high repeatability of results,
- LCD display of actual parameters: speed, time and relative centrifugal force. Optimization of acceleration and deceleration force.
- "Short" spinning function.
- Inspection glass.
- Unbalance protection system.
- It is not possible to open the lid, during spinning.
- It is not possible to start, when the lid is open.

Product according to national and international safety regulations EN-61010-1 and EN-61010-2-020 for laboratory equipments (CE).

Ref	Dimensions W x H x D	Weight	Consumption W
<b>1460</b>	355x270x435mm.	13Kg	170
Capacity	Max. R.p.m	Max. RCF	
500 ml.	5000	3300 xg.	
Maintenance-free	Timer	Voltage	Brussless
Yes	1-99 min ∞	130/250v. 50/60Hz.	Yes
Display	Rotor identification	Refrigerant	Temperature range
LCD	NO	NO	---

Rotor		Accessories information				
 Ref. 1460-1487C/1237C	8 y 12x15ml.	Accessories	Ø x Long	Tubes	Ref.	
	Angle 45/30 Max. 5000 R.p.m. Max. 3300 xg	Round carriers 5ml.	Ø13x75mm.	8 / 12	<b>4088</b>	
		Round carriers 7ml.	Ø13x100mm	8 / 12	<b>4082</b>	
Rotor		Accessories information				
 Ref. 1460-1486C	24x15ml.	Accessories	Ø x Long	Tubes	Ref.	
	Angle 30 Max. 4000 R.p.m. Max. 2320 xg	Round carriers 5ml.	Ø13x75mm.	24	<b>4088</b>	
		Round carriers 7ml.	Ø13x100mm	24	<b>4082</b>	
Rotor		Accessories information				
 Ref. 1460-1257C	4x25/30ml.	Accessories	Ø x Long	Tubes	Ref.	
	Angle 45 Max. 5000 R.p.m. Max. 3300 xg		Ø25x100mm	4		
		Round carriers 5 / 7ml.			4	<b>4255</b>
		Round carriers 15 /10ml.			4	<b>4256</b>
Rotor		Accessories information				
 Ref. 1460-1329C/1332C	10 y 12x25/30ml.	Accessories	Ø x Long	Tubes	Ref.	
	Angle 40/30 Max. 4000 R.p.m. Max. 2320 xg		Ø25x100mm	10/12		
		Round carriers 5 / 7ml.			10/12	<b>4255</b>
		Round carriers 15 /10ml.			10/12	<b>4256</b>
Rotor		Accessories information				
 Ref. 1460-1333C	10x50ml. Fal.	Accessories	Ø x Long	Tubes	Ref.	
	Angle 30 Max. 4000 R.p.m. Max. 2320 xg	Buckets	Ø35x120mm	10	<b>3276</b>	
		Round carriers 30 / 25ml.	Ø25x120mm	10	<b>4248</b>	

Rotor		Accesories information			
	4x100ml.	Accesories	Ø x Long	Tubes	Ref.
	Swing out Max. 4000 R.p.m. Max. 2469 xg	Hangers with buckets for 15ml. Fal	Ø17x120mm	8	<b>3195C</b>
		Buckets for 100ml.	Ø44.5mm	4	<b>3182</b>
		Buckets W/C for 100ml.	Ø44.5mm	4	<b>3184C</b>
		Round carriers 50ml. Fal.	Ø30x100mm	4	<b>4189</b>
		Round carriers 50ml.	Ø35x100mm	4	<b>4192</b>
		Round carriers 25 / 30ml.	Ø25x100mm	4	<b>4190</b>
		Round carriers 15 / 10ml.	Ø16x100mm	16	<b>4187</b>
		Round carriers 5 / 7ml.	Ø13x75/100mm	16	<b>4186</b>
		Cups alum. Secur. Biolog.		4	<b>7185</b>
Rubber under tube 100/50/30/25ml.		4	<b>4188</b>		
Ref. 1460-2183					
Rotor		Accesories information			
	8x15ml.	Accesories	Ø x Long	Tubes	Ref.
	Swing out Max. 4000 R.p.m. Max. 1967 xg				
		Round carriers 5ml.	Ø13x75mm.	8	<b>4088</b>
		Round carriers 7ml.	Ø13x100mm	8	<b>4082</b>
Ref. 1460-2193C					
Rotor		Accesories information			
	4x70ml	Accesories	Ø x Long	Tubes	Ref.
	Swing out Max. 4000 R.p.m. Max. 2320 xg	Buckets	Ø70mm	4	<b>3233</b>
		Hangers for 5ml.	Ø13x81mm	48	<b>3170</b>
		Round carriers 5ml.	Ø13x81mm	32	<b>4169</b>
		Round carriers 5ml. Sin tapón	Ø12x75mm	48	<b>4235</b>
		Round carriers 5ml.	Ø13x75mm	28	<b>4236</b>
		Round carriers 10 / 7ml.	Ø17x75mm	28	<b>4234</b>
Ref. 1460-2232					

Rotor		Accessories information			
 Ref. 1460-2485	4x50ml.	Accessories	Ø x Long	Tubes	Ref.
	Swing out Max. 4000 R.p.m. Max. 2320 xg	Hangers for 5ml.	Ø13x75mm	20	<b>3245</b>
		Hangers with buckets for 15ml.	Ø16x100/110mm	16	<b>3246C</b>
		Hangers with buckets for 15ml. Fal	Ø17x120mm	8	<b>3247C</b>
		Buckets W/C 50ml. Fal.	Ø35x100mm	4	<b>3249C</b>
		Buckets 50ml. Fal.	Ø35x100mm	4	<b>3250</b>
	Round carriers 30 / 25ml.	Ø25x100mm	4	<b>4248</b>	
Rotor		Accessories information			
 Ref. 1460-2331C	4x10ml.	Accessories	Ø x Long	Tubes	Ref.
	Ø17x85mm				
	Swing out Max. 2500 R.p.m. Max. 750 xg	Round carriers 5ml.	Ø13x75mm.	4	<b>4088</b>
		Round carriers 7ml.	Ø13x100mm	4	<b>4082</b>
Rotor		Accessories information			
 Ref. 1460-2270C	12 Cards	Accessories	Ø x Long	Tubes	Ref.
	Column test Swing out Max. 1000 R.p.m. Max. 85 xg				