



Ref. DA 1040 – DA 1048

GENERAL CHARACTERISTICS

- Designed for the distilled water production from local supplying. Easy cleaning and security systems.
- Double wall, metal construction: exterior painted in epoxy and interior made of stainless steel.
- Sealed with silicon.
- Control panel with general interrupter and temperature control.
- Security system with hydraulic temperature control.
- Liquids level sensor with automatic cycle stopping in case of lack of water
- Water entry connection adaptable to the entry pipe.
- Exit connector adaptable to a container.

Manufactured according to EC directives

Ref.	Capacity	Dimensions W x H x D	Weight.	Consumption	Cooling water
DA1040	4 L/h.	250 x 540 x 220mm.	12 Kg.	3000 W.	30 L.
DA1048	8 L/h.	260 x 610 x 260mm.	14 Kg.	6000 W.	60 L.