



Ref. 1150

GENERAL CHARACTERISTICS

- Designed for continuous and smooth rotation of blood bags.
- It consists of a stainless steel basket with 36 compartments connected to a motor that rotates the basket continuously at very slow speed.
- Metal housing painted in epoxy.
- Rotation capacity: 36 bags of 500 ml.
- AISI 304 stainless steel basket.
- Illuminated ON indicator light.
- It incorporates a circuit breaker located at the back of the motor housing.

Manufactured according to CE standards.

Ref.	Rotation speed (rpm)	Consumption	Dimensions (mm) W x H x D	Weight
1150	Fixed: 8	40 W	700 x 250 x 300	6 Kg.