



Ref. 1440



Ref. 1420

### GENERAL CHARACTERISTICS

Shaker with vibrating movement for shaking four 85 x 125 mm microtiter plates or any other kind of plates, given that trays can be interchanged (ref.1440). There is also another version that allows adapting flasks of many different sizes (ref.1420).

- Electronic speed adjustment.
- Great stability at maximum shaking.
- 60 minutes timer with a continuous switch function.

Manufactured according to CE directives

Ref.	Vibration orbit	Speed (rpm)	Consumption	Max. load	Useful platform	Dimensions W x H x D	Weight
1440	3 mm.	Adjustable from 100-1000	20 W	2,5 Kg	220 x 220	230 x 115 x 295	6 Kg
1420	6 mm.	Adjustable from 100-500	20 W	2,5 Kg	220 x 220	230 x 115 x 295	6 Kg

Ref.	Accesorios
1450-1	50ml. Flask adaptor
1450-2	100ml. Flask adaptor
1450-3	250ml. Flask adaptor
1450-4	500ml. Flask adaptor
1450-5	1000ml. Flask adaptor

