



Ref. 1730

GENERAL CHARACTERISTICS

- Electronic regulator that allows variation of the power of the heating elements by means of semiconductors.
- Electronic regulation of the voltage that allows variation of the heating power between 10% and 100% of this voltage.
- Incorporates a suppression filter for interferences.
- Incorporates a power outlet for fixing a contact thermometer.

Manufactured according to EC directives

Ref.	Dimensions W x H x D	Maximun charge permitted	Weight
1730	200 x 100 x 190mm.	2500W – 230V.	1,5 Kg.