

Ref. 1110



Ref. 1100



GENERAL CHARACTERISTICS

- Head stock for mixing test tubes of up to Ø 20 mm.
- It has two switchable working modes:
 - Continuous mode for uninterrupted mixing.
 - Manual mode in which mixing is activated automatically by slight pressure from the tube on the head stock.
- Metal housing painted in epoxy.
- Electronically regulated speed control.

Manufactured according to CE directives

Ref.	Mixing speed (rpm)	Consumption (W)	Dimensions (mm) W x H x D	Weight (Kg)	Eccentricity
1100	between 50 and 2600	48	140 x 100 x 155	4	6mm.
1110	Fixed: 2600	48	140 x 100 x 155	4	6mm.



Ref.	Accessories
1110-1	Tube rack for 16 test tubes (ømax 15mm)
1110-2	Tube rack for 25 test tubes (ømax 12mm)
1110-3	Head stock for 6 Eppendorf tubes
1110-4	90-mm platform for containers of up to 250 ml

