

Ref. 1321



Ref. 1320



GENERAL CHARACTERISTICS

- Laboratory stirrer for stirring jobs from 10 to 20 litres (relative to water).
- Metal housing painted in epoxy.
- Universal clamp for rapid coupling of rods of up to Ø10 mm.
- Fastening system located at the back of the equipment.
- Electronic speed control variable from 200 to 2000 r.p.m.
- A standard, Ø 8 mm, 300 mm long stirring rod and a 80 x 20 mm blade is supplied with the AGV-10.

Manufactured according to CE directives.

Ref.	RPM indicator	Power	Maxtorque(Ncm)	Weight	Stirring capacity (l)	Dimensions W x H x D
1320	Not included	50 W	75	4,5 Kg.	Up to 20	90 x 200 x 180mm
1321	Digital with reading precision of 2%	50 W	75	4,5 Kg.	Up to 20	90 x 200 x 180mm



Ref.	Accessories
1322	Support foot for Ø 15 x 750mm bar and chuck
1323	Standard rod
1324	Crossed blade rod
1325	Paddle-shaped rod
1326	Drop rod