

Agitador para Bolsas de Sangre SERIE GR-2

Blood Bag Rotator GR-2 SERIES



Ref. 1150

CARACTERÍSTICAS GENERALES

- Diseñado para una agitación continua y suave de bolsas de sangre.
- Consta de un cestillo de acero inoxidable con 36 departamentos. Indicado para bancos de sangre que permite, dado su constante movimiento, impedir que la sangre se coagule, permitiendo una uniformidad en las muestras contenidas en las bolsas.
- Mueble metálico pintado en epoxi.
- Capacidad de agitación: 36 bolsas de 500 ml.
- Cestillo de acero inoxidable AISI 304.
- Interruptor luminoso de puesta en marcha.
- Incorpora fusible de protección situado en la parte posterior de la caja motor.

Fabricado según normas CE.

GENERAL CHARACTERISTICS

- Designed for continuous and smooth rotation of blood bags.
- It consists of a stainless steel basket with 36 compartments.
- Metal housing painted in epoxy.
- Rotation capacity: 36 bags of 500 ml.
- AISI 304 stainless steel basket.
- Illuminated ON indicator light.
- It incorporates a fuse protection located at the back of the motor housing.

Manufactured according to CE standards.

Ref.	Velocidad de rotación (rpm) Rotation speed (rpm)	Potencia (W) Consumption (W)	Dimensiones (An x Al x F, mm) Dimensions (W x H x D, mm)	Peso (Kg) Weight (Kg)
1150	8	40	700 x 250 x 300	6