

Agitador Miniorbital SERIE MO-200

Miniorbital shaker MO-200 SERIES

CARACTERÍSTICAS GENERALES

Agitador con movimiento vibratorio para agitar 4 placas microtiter de 85 x 125 mm ó cualquier otro tipo de placas debido a la posibilidad de intercambiar bandejas (ref.1440). También existe otra versión con la posibilidad de adaptar un sistema de fijación universal para diferentes recipientes (ref.1420).

- Regulación electrónica de la velocidad.
- Gran estabilidad a máxima agitación.
- Reloj temporizador; seleccionable hasta 60 minutos, con posición para funcionamiento continuo.

Fabricado bajo directivas CE.

GENERAL CHARACTERISTICS

Shaker with vibrating movement for shaking four 85 x 125 mm microtiter plates or any other kind of plates, given that trays can be interchanged (ref.1440). There is also another version that allows adapting flasks of many different sizes (ref.1420).

- Electronic speed adjustment.
- Great stability at maximum shaking.
- 60 minutes timer with a continuous switch function.

Manufactured according to CE directives.



Ref. 1440



Ref. 1420

Ref.	Órbita de vibración (mm) Vibration orbit (mm)	Velocidad (rpm) Speed (rpm)	Consumo (W) Consumption (W)	Carga máx. (Kg) Max. load (Kg)	Plataforma útil Useful platform	Dimensiones (An x Al x F, mm) Dimensions (W x H x D, mm)	Peso (Kg) Weight (Kg)
1440	3	Regulable de 100-1000 Adjustable from 10-1000	20	2,5	220 x 220	230 x 115 x 295	6
1420	6	Regulable de 100-500 Adjustable from 100-500	20	2,5	220 x 220	230 x 115 x 295	6

Ref.	Accesorios Accessories
1450-1	Adaptador para matraces de 25/50 ml. / 25/50 ml. Flask adaptor
1450-2	Adaptador para matraces de 100 ml. / 100 ml. Flask adaptor
1450-3	Adaptador para matraces de 250 ml. / 250 ml. Flask adaptor
1450-4	Adaptador para matraces de 500 ml. / 500 ml. Flask adaptor
1450-5	Adaptador para matraces de 1000 ml. / 1000 ml. Flask adaptor