



Ref. 1730

CARACTERISTICAS GENERALES

- Regulador electrónico que permite variar la potencia de los elementos calefactores por medio de semiconductores.
- Regulación electrónica de la tensión que permite variar la potencia de calefacción entre un 10% y un 100% de esta.
- Incluye un filtro supresor de interferencias.
- Incorpora una toma para colocar termómetro de contacto.
- Fabricado bajo directivas CE
- Tensión: 230v/50H

GENERAL CHARACTERISTICS

- Electronic regulator that allows variation of the power of the heating elements by means of semiconductors.
- Electronic regulation of the voltage that allows variation of the heating power between 10% and 100% of this voltage.
- Incorporates a suppression filter for interferences.
- Incorporates a power outlet for fixing a contact thermometer.
- Manufactured according to CE directives
- Tension: 230V/50Hz

Ref.	Carga máxima permitida Maximum load permitted	Dimensiones ø ext. x al Dimensions ext. ø x H (mm)	Peso Weight (Kg)
1730	2500W - 230V	200 x 100 x 190	1,5