

Agitador Orbital y de Vaivén SERIE AO-400

Linear and Orbital shaker AO-400 SERIES



Ref. 1400

CARACTERÍSTICAS GENERALES

- Movimiento orbital o de vaivén modificable por el usuario.
- Chasis robusto que permite altas velocidades de agitación con gran estabilidad.
- Mueble metálico pintado en epoxi.
- La bandeja de agitación está diseñada para acoplar distintos recipientes de laboratorio y está provista de cuatro barras metálicas acolchadas para fijar de varias formas dichos recipientes.
- Entre sus especificaciones técnicas destaca:
 - Velocidad de agitación; variable de 40 a 300 r.p.m. mediante regulación electrónica.
 - Reloj temporizador; seleccionable hasta 60 minutos, con posición para funcionamiento continuo.
 - Bandeja estándar para agitación de 450 x 350 mm. (ancho x fondo).
 - Cuentarrevoluciones digital con una precisión del $\pm 2\%$.

Fabricados según normas CE.



Ref. 1400

GENERAL CHARACTERISTICS

- Equipment designed for reproducible shaking of flasks, Erlenmeyer flasks, bottles, test tubes, decantation funnels, "Petri" capsules, etc.
- Orbital o linear movement.
- A strong chassis that allows high shaking speeds with great stability.
- Metal housing painted in epoxy.
- The shaking platform is designed for coupling different laboratory containers, and it is equipped with four, cushioned metal bars for fastening containers in various ways.
- The following are highlights of its technical specifications:
 - Shaking speed: variable from 40 to 300 r.p.m. by electronic regulation.
 - Timer: selectable up to 60 minutes, with a continuous operation position.
 - Standard, 450 x 350 mm platform (width x depth) for shaking.
 - Digital R.P.M. indicator with a precision of $\pm 2\%$.

Manufactured according to CE standards.

Ref.	Capacidad de agitación (Kg) Shaking capacity (Kg)	Consumo (W) Consumption (W)	Excentricidad Eccentricity	Dimensiones sin bandeja (An x Al x F) Dimensions w/o tray (W x H x D)	Peso (Kg) Weight (Kg)
1400	Hasta 15 Up to 15	90	15 mm.	435 x 140 x 365 mm.	18

Agitador Orbital y de vaivén SERIE AO-400

Linear and Orbital shaker AO-400 SERIES



Ref. 1450



Ref. 1450-5



Ref. 1450-4



Ref. 1450-3



Ref. 1450-2



Ref. 1450-1



Ref.	Accesorios Accessories
1450	Plataforma de 350 x 450 mm para matraces Erlenmeyer 400x500 mm platform for Erlenmeyer flask
1450-1	Adaptador para matraces de 25/50 ml. (42 adaptadores en la plataforma) 25/50 ml. Flask adaptor (capacity for up 42 adaptors on each platform)
1450-2	Adaptador para matraces de 100 ml. (20 adaptadores en la plataforma) 100 ml. Flask adaptor (capacity for up 20 adaptors on each platform)
1450-3	Adaptador para matraces de 250 ml. (12 adaptadores en la plataforma) 250 ml. Flask adaptor (capacity for up 12 adaptors on each platform)
1450-4	Adaptador para matraces de 500 ml. (9 adaptadores en la plataforma) 500 ml. Flask adaptor (capacity for up 9 adaptors on each platform)
1450-5	Adaptador para matraces de 1000 ml. (6 adaptadores en la plataforma) 1000 ml. Flask adaptor (capacity for up 6 adaptors on each platform)
	Bandejas especiales bajo demanda Special platform upon request