



Ref. 1440



Ref. 1420

CARACTÉRISTIQUES GÉNÉRALES

Agitateur à mouvement rotatif. Son plateau accepte tout type de plaque (réf : 1440) comme par exemple, 4 plaques de 85 x 125 mm pour la micro-titration et il peut être équipé pour recevoir d'autres types de flacons de différentes tailles (réf : 1420)

- Régulation électronique de la vitesse.
- Grande stabilité à la vitesse maximale.
- Minuterie 0 – 60 min et mode de fonctionnement continu
- Fabriqué en conformité avec les normes CE
- Alimentation : 230 V – 50 Hz

GENERAL CHARACTERISTICS

Shaker with vibrating movement for shaking four, 85 x 125-mm microtiter plates or any other kind of plates, given that trays can be interchanged (ref.1440). There is also another version that allows adapting flasks of many different sizes (ref.1420)

- Electronic speed adjustment.
- Great stability at maximum shaking.
- 60 minutes timer with a continuous switch function
- Manufactured according to CE directives
- Voltage: 230V/50Hz

Ref.	Amplitude de l'orbite (mm) Vibration orbit (mm)	Vitesse (rpm) Speed (rpm)	Puissance (W) Consumption (W)	Charge máx. admissible (Kg) Max. allowable load (Kg)	Plateau utile (mm) Useful platform (mm)	Dimensions (mm) L x H x P Dimensions (mm) W x H x D	Poids (Kg) Weight (Kg)
1440	3	Réglable de 100-1000 Adjustable from 100-1000	20	2,5	220 x 220	230 x 115 x 295	6
1420	6	Réglable de 100-500 Adjustable from 100-500	20	2,5	220 x 220	230 x 115 x 295	6

Ref.	Accesorios
1450-1	Adaptateur pour flacon 50 ml.
1450-2	Adaptateur pour flacon 100 ml.
1450-3	Adaptateur pour flacon 250 ml.
1450-4	Adaptateur pour flacon 500 ml.
1450-5	Adaptateur pour flacon 1000 ml.

Ref.	Accessories
1450-1	50ml. Flask adaptor
1450-2	100ml. Flask adaptor
1450-3	250ml. Flask adaptor
1450-4	500ml. Flask adaptor
1450-5	1000ml. Flask adaptor



ref. 1450-5

ref. 1450-4

ref. 1450-3

ref. 1450-2

ref. 1450-1