



[www.bunsen.es](http://www.bunsen.es)



# ANALOG HEATING MANTLE

Ref. 1703

## GENERAL CHARACTERISTICS

- \* Electronic power regulator
- \* Illuminated operating lights.
- \* ABS case.
- \* Excellent thermal insulation with fibroceramic.

Manufactured according to EC directives.

| Ref  | Capacity | Dimensions Ø x H | Useful measures Ø x H | Max. temp | Power (W) |
|------|----------|------------------|-----------------------|-----------|-----------|
| 1703 | 1000ml.  | Ø260 x 165mm.    | Ø145 x 85mm.          | 380°C     | 530       |



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fa. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

