



www.bunsen.es



DIGITAL HEATING MANTLE

Ref. 1711

GENERAL CHARACTERISTICS

- * Temperature: Digital controller with PID control.
- * Digital display with real and programmed temperature.
- * External probe for temperature control in sample.
 - * ABS Case.
- * Excellent thermal insulation with fibroceraamic.

Manufactured according to EC directives.

Ref	Capacity	Dimensions Ø x H	Useful measures Ø x H	Max. temp	Power (W)
1711	250ml.	Ø230 x 155mm.	Ø95 x 60mm.	380°C	210



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

