



www.bunsen.es



DIGITAL HEATING MANTLE

Ref. 1713

GENERAL CHARACTERISTICS

- * Temperature: Digital controller with PID control.
- * Digital display with real and programmed temperature.
- * External probe for temperature control in sample.
 - * ABS Case.
- * Excellent thermal insulation with fibroc ceramic.

Manufactured according to EC directives.

Ref	Capacity	Dimensions Ø x H	Useful measures Ø x H	Max. temp	Power (W)
1713	1000ml.	Ø260 x 165mm.	Ø145 x 85mm.	380°C	530



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

