



www.bunsen.es



CIRCULAR HEATING PLATE

Ref. 1502

CARACTERÍSTICAS GENERALES

- * Circular heating plate of solid steel Ø 180mm. Power 1500 W.
 - * 7-power switch.
 - * External cabinet 265 x 265 mm.
- * Enamelled box with high resistance to heat and chemicals.
 - * Pilot operating indicator.

Manufactured according to EC directives.



Ref	Dimensions W x H x D	Surface	Max. temp	Power (W)	Weight (Kg)
1502	265 x 125 x 265mm.	Ø 150 mm.	400°C	1500	2



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

