




[www.bunsen.es](http://www.bunsen.es)



# PLATELETS SHAKER

Ref. 1400-BS

## GENERAL CHARACTERISTICS

Equipment designed for continuous shaking of platelets or similar bags.

- \* Metal housing painted in epoxy.
- \* Linear movement, 12mm lenght.
- \* Shaking speed, fixed 60 rpm.
- \* Special rack to shake up to 18 bags (stainless steel rack).
- \* Self-oiling sliding system.

Manufactured according to CE standards.

Ref	Stirring capacity	Consumption (W)	Dimensions (W x H x D, mm)	Weight (Kg)	Rack dimensions (mm)
1400-BS	18 bolsas	120	505 x 425 x 430	20	385 x 285 x 355



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fax. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

