



[www.bunsen.es](http://www.bunsen.es)



# TUBE ROTATOR

Ref. 1130

## GENERAL CHARACTERISTICS

This rotator is especially designed with the following characteristics for slowly and continuously rotating test tubes.

- \* Injected aluminium housing painted in epoxy.
- \* Ø 280 mm test tube adapter disk. **NO INCLUDING.**
- \* Possibility of varying the angle of rotation axis from 0 to 90 degrees.
  - \* Electronically regulated speed control.
- \* Possibility of shaking test tubes from Ø 6 mm to Ø 38 mm.

Manufactured according to CE standards.



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf.+34 916 113 584 - Fax. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)



Ref	Rotation speed (rpm)	Consumption (W)	Dimensions (W x H x D, mm)	Weight (Kg)
<b>1130</b>	Variable between 5 y 25	12	180 x 300 x 220	4

Ref	Accessories
<b>1130-1</b>	Disks for 24 tubes from 10 to 12 mm
<b>1130-2</b>	Disks for 20 tubes from 14 to 16 mm
<b>1130-3</b>	Disks for 20 tubes from 16 to 20 mm
<b>1130-4</b>	Disks for 16 tubes from 20 to 24 mm
<b>1130-5</b>	Disks for 16 tubes from 26 to 32 mm
<b>1130-6</b>	Disks for 16 tubes from up to 40 mm
	Disks, test tubes and vessels with other dimensions available on request



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf.+34 916 113 584 - Fax. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

