



www.bunsen.es



ROD STIRRER

Ref. 1320

GENERAL CHARACTERISTICS

- * Shaking capacity, 20L. (related to water).
- * Metal housing painted in epoxy.
- * Universal clamp for rapid coupling of rods of up to Ø 10 mm.
- * Double nut fastening system located on the back of the equipment.
- * Speed; 200 ÷ 2000 r.p.m. adjustable by means of electronic regulation.
- * Supplied with standard shaking rod of 300mm of length, 8mmØ and an 80x20mm. blade, of AISI 304 stainless steel.
- * Maximum viscosity: 50000 mPas.

Manufactured according to CE directives.



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es



Ref	Max. torque (Ncm)	Power (W)	Dimensions (W x H x D, mm)	Weight (Kg)	Stirring capacity
1320	27	90	105 x 160 x 220	4.5	Up to 20L.

Ref	Accessories
1322	Support foot for Ø 20 x 750 mm bar and chuck
1323	Standard rod
1324	Crossed blade rod
1325	Paddle/shaped rod
1326	Drop rod



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es

