



www.bunsen.es



DIGITAL ROD STIRRER

Ref. 1301

GENERAL CHARACTERISTICS

- * Metal housing painted in epoxy.
- * Universal clamp for rapid coupling of rods of up to Ø10mm.
 - * Microprocessor controlled.
- * Display of speed and time settings using backlit LCD display.
 - * Digital rev counter with accuracy of $\pm 2\%$.
- * Timed mode up 9h. 59'59". (Factory-changeable).
 - * Fixing system located on the back of the unit.
 - * Adjustable speed 40 ÷ 2000 rpm.
- * Stable speed throughout the range, regardless of load (Max. load 10L. Related to water).
 - * Supplied with standard shaking rod of 250 mm of length, 8 mm Ø and 50x10 mm blade, of AISI 304 stainless steel.
 - * Maximum Viscosity: 1000 mPas.

Manufactured according to CE directives.



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es



Ref	Max. torque (Ncm)	Power (W)	Dimensions (W x H x D, mm)	Weight (Kg)	Stirring capacity
1301	5.3	20	120 x 145 x 160	3.5	Up to 10L.

Ref	Accessories
1312	Support foot for Ø 15 x 750 mm bar and chuck
1313	Standard rod
1314	Crossed blade rod
1315	Paddle/shaped rod



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es

