



[www.bunsen.es](http://www.bunsen.es)



# BLOOD BAG ROTATOR

Ref. 1150

## GENERAL CHARACTERISTICS

- \* Designed for continuous and smooth rotation of blood bags.
- \* It consists of a stainless steel basket with 36 compartments.
  - \* Metal housing painted in epoxy.
  - \* Rotation capacity: 36 bags of 500 ml.
  - \* AISI 304 stainless steel basket.
  - \* Illuminated ON indicator light.
- \* It incorporates a fuse protection located at the back of the motor housing.

Manufactured according to CE standards.

Ref	Rotation speed (rpm)	Consumption (W)	Dimensions (W x H x D, mm)	Weight (Kg)
1150	8 Rpm	40	700 x 250 x 300	6



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fax. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

