



[www.bunsen.es](http://www.bunsen.es)



## DIGITAL BATH 12L. 200°C

Ref. 1614

### GENERAL CHARACTERISTICS

- \* Double wall, metal construction: exterior painted in epoxy with antioxidant treatment and interior basin made of 18/10 stainless steel.
  - \* Heating by shielded heating element made of AISI 316 stainless steel.
- \* Equipped with a safety thermostat adjustable that has a manual reset in the event of over-temperature according to DIN 12877, class 2.
- \* Digital Temperature control by digital thermostat controlled by microcontroller from environment + 5° to 200°C. Resolution, 0.1°C.
- \* Temperature regulation by means of a digital thermostat with auto-adjustable PID control (AUTOTUNIG), in order to obtain very low working temperature inertia.
  - \* Temperature display with double digital display of 4 digits.
- \* Control panel includes:
  - \* Illuminated ON switch.
  - \* Indicator light showing heating element operation.
  - \* Safety thermostat indicator light.
  - \* Display control.

Manufactured according to CE directives.



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fax. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)



Ref	Capacity	Probe	Power	Useful measures W x H x D	Ext. measures W x H x D	Max. temp.	Precision °C
<b>1614</b>	12L.	PT-100	1500W.	320 x 150 x 290mm.	390 x 200 x 440mm.	200°C	±0.5%

Ref	Accessories
<b>1631</b>	Stainless steel cover with 4 openings and reduction disks
<b>1641</b>	Top cover
<b>1651</b>	Element cover
<b>1661</b>	Blank cover
<b>1663</b>	Drain tap
<b>1680</b>	Hinge for top cover ref. 1641
<b>CHB 061</b>	Racks for baby bottles (12 baby bottles)
<b>1600-10-1</b>	Racks (265x120 mm) for 98 test tubes, Ø 10-12 mm
<b>1600-10-2</b>	Racks (265x120 mm) for 72 test tubes, Ø 14-16 mm
<b>1600-10-3</b>	Racks (265x120 mm) for 50 test tubes, Ø 18-20 mm
<b>1600-3</b>	Liquid level indicator (factory installed)
	Different platform sizes upon client's request



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN  
Tlf. +34 916 113 584 - Fax. +34 916 128 254  
[info@bunsen.es](mailto:info@bunsen.es) - [www.bunsen.es](http://www.bunsen.es)

