




www.bunsen.es

DRYING OVEN



115L. DRYING OVEN WITH FORCED AIR CONVECTION Ref. 8311

GENERAL CHARACTERISTICS

- * Interior and exterior construction in stainless steel AISI 304.
- * Solid door, fitted with a key. Excellent thermal insulation.
- * Protection against overtemperature class 2.0 according to DIN 12880.
 - * 4,3", clear, full colour touch screen.
 - * USB and LAN ports for data transfer
- * Multi-segment time and temperature programs. Loop function up to 255 times or endless.
 - * Adjustable start delay feature, adjustable ramps.
- * Recording of min, average and max temperature value for each segment overview of set and current parameters while operating.
 - * Operating in temperature or time priority mode.
- * Fan speed control in the range of 0 ... 100% or 10% ... 100% (depending on the chamber capacity).
 - * Temperature sensor fail alarm.
 - * Possibility of temperature calibration by the User.
- * Power failure control system (program continued after restoring power).
 - * Visual and sound alarm.
 - * Quick change of program parameters.



BUNSEN SA. 30 Rio Alberche St - 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es













* Alarm Bar.

* Automatic fan shut-down after completing the program. Automatic fan air-flap control.

* Shelves standard included 2, possibility of including 7. Max shelf workload 25Kg.

Manufactured according to CE directives.

Ref.	Capacity	Dimensions W x H x D	Useful measures (W x H x D)	Weight	Temperature	Nominal power
8311	115L.	650 x 850 x 710mm.	460 x 540 x 450mm.	65 Kg.	Up to 300°C	2500

Ref	Accessories	
8300-02		Stainless steel wire shelf INOX
8300-03		Perforated shelf
8300-04		Reinforced shelf
8300-05		Reinforced version
8300-06	1, 2, 3, ... 	Door openings counter
8300-07		Table with wheels
8300-08		Additional Pt 100 temperature sensor
8300-09		HEPA-fresh air filter
8300-10		Wheels
8300-11		Non-standard access port for external sensor



BUNSEN SA. 30 Rio Alberche St - 28970 Humanes de Madrid, Madrid, SPAIN

Tlf. +34 916 113 584 - Fax. +34 916 128 254

info@bunsen.es - www.bunsen.es

