



www.bunsen.es

Low SPEED



GEN PLUS CENTRIFUGE

Ref. 1469

GENERAL CHARACTERISTICS

- * Maintenance – free, brush less motor.
- * Spinning in “short” mode.
- * Controlled by microprocessor with LCD screen that permits displaying all parameters.
 - * Timer: 1 ÷ 99 min. or ∞.
- * Rotor chamber made from stainless steel.
 - * Inspection glass.
- * Safety features:
 - * it is not possible to open the lid during spinning.
 - * it is not possible to start when the lid is open.
 - * Unbalance sensor.
 - * Emergency lid lock release.

Manufactured according to National and International standards EN-61010-1 and EN-61010-2-020 on safety and EC marking.



BUNSEN SA. C/ Rio Alberche, 30 - 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es




Ref	Dimensions W x H x D	Weight	Power	Max. Capacity	Max. Speed	Max. RCF
1469	380 x 430 x 540mm.	43 Kg.	320W.	1000mL.	100 ÷ 4500 Rpm.	3305 (xg)
Timer (min)	Voltage		Display	Rotor ID	Refrigeration	Temperature Range
0-99	115 v/ 230 v -- 50/60Hz		LCD	Si	No	-----

Ref	Accessories	
1469-1210C		Angular 24 x 15 mL. / Max: 4500 Rpm. / RCF Max: 3237 xg / 30°
1469-1465C		Angular 30 x 15 mL. / Max: 4500 Rpm. / RCF Max: 3011 xg / 30°
1469-1215C		Angular 32 x 15 mL. / Max: 4500 Rpm. / RCF Max: 3283 xg / 45°
1469-1743C		Angular 12 x 30 mL. / Max: 4500 Rpm. / RCF Max: 2716 xg / 30°
1469-3275C		Angular 10 x 50mL. / Max: 4500 Rpm. / RCF Max: 2416 xg / 30° Biological security.
1469-3276		Angular 10 x 50mL. / Max: 4500 Rpm. / RCF Max: 2416 xg / 30°
1469-2436a		Swing out 4 x 200 mL. / Max: 4500 Rpm. / RCF Max: 3305 xg.
1469-2436c		Swing out 4 x 200 mL. / Max: 4500 Rpm. / RCF Max: 3305 xg. Biological security.



BUNSEN SA. C/ Rio Alberche, 30 - 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es



Ref	Accesorios	
1469-2177		Swing out 4 x 250ml. / Max: 4500 Rpm. / RCF Max: 4150 xg.
1469-2177C		Swing out 4 x 250ml. / Max: 4500 Rpm. / RCF Max: 4150 xg. Biological security.
1469-2451C		Swing out 6 x P. microtiter. / Max: 3000 Rpm. / RCF Max: 1036 xg.

Más centrifugas



Other centrifuges



BUNSEN SA. C/ Rio Alberche, 30 - 28970 Humanes de Madrid, Madrid, SPAIN

Tlf. +34 916 113 584 - Fax. +34 916 128 254

info@bunsen.es - www.bunsen.es

