



www.bunsen.es

High SPEED



HEMO PLUS CENTRIFUGE

Ref. 1461

GENERAL CHARACTERISTICS

- * Controlled by microprocessor with large graphic LCD - blanking the screen after a period of inactivity.
 - * Two screen modes (simplified and normal) with 10 menu languages.
 - * 100 user programs and 10 acceleration/deceleration curves.
 - * Speed or RCF adjustment.
 - * Manual centrifugal radius adjustment (with automatic RCF correction).
- * Time calculation from start key pressing or from reaching preselected speed. Ascending or descending time counting.
 - * Changing of parameters during centrifuging.
 - * Visual and sound alerts of different notifications.
 - * Automatic rotor recognition.
 - * Selected functions blocking, password protection.
- * Continuous operation mode HOLD and SHORT time operation mode.
 - * Automatic lid opening.



BUNSEN SA. C/ Rio Alberche, 30 - 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es



- * Selected functions blocking, password protection.
- * Stainless steel rotor chamber. Unbalance sensor.
- * Lid locking during rotor running. Emergency lid lock release.
- * Start blocking at opened lid.

Manufactured according to National and International standards EN-61010-1 and EN-61010-2-020 on safety and EC marking.


Ref	Dimensions W x H x D	Weight	Power	Max. Capacity	Max. Speed	Max. RCF
1461	367 x 300 x 455	22 Kg.	220W.	500mL.	90 ÷ 18000 Rpm.	24270 (xg)
Timer (min)		Voltage	Display	Rotor ID	Refrigeration	Temperature Range
1s ÷ 99h 59min 59s, ∞		115v 50/60Hz/230v 50/60Hz	LCD	Yes	No	-----


Ref	Accessories	
1461-1199		Angular 12 x 1.5/2.2 mL. / Max: 18000 Rpm. / RCF Max: 24088 xg 45°. Hermetically sealed lid.
1461-1461		Angular 24 x 1.5/2.2 mL. / Max: 14000 Rpm. / RCF Max: 18187 xg 45°. Hermetically sealed lid.
1461-1462		Angular 36 x 1.5/2.2 mL. / Max: 14000 Rpm. / RCF Max: 18187 xg 45°. Hermetically sealed lid.
1461-1716		Angular (4x8) x 0.2 mL. / Max: 14000 Rpm. / RCF Max: 15339 xg / 45°. Hermetically sealed lid. PCR strips.
1461-1216		Angular 12 x 5 mL. / Max: 14000 Rpm. / RCF Max: 19064 xg / 30°. Hermetically sealed lid.
1461-1715		Angular 10 x 10 mL. / Max: 14000 Rpm. / RCF Max: 15558 xg / 30°. Hermetically sealed lid.



BUNSEN SA. C/ Rio Alberche, 30 - 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es



Ref	Accessories	
1461-1741C		Angular 8 x 10/15 mL. / Max: 6000 Rpm. / RCF Max: 4226 xg / 30°.
1461-1217C		Angular 10 x 10/15 mL. / Max: 6000 Rpm. / RCF Max: 4226 xg / 30°.
1461-1740C		Angular 12 x 10/15 mL. / Max: 6000 Rpm. / RCF Max: 4226 xg / 30°.
1461-1745C		Angular 24 x 10/15 mL. / Max: 5000 Rpm. / RCF Max: 3354 xg / 30°.
1461-1501C		Angular 30 x 10/15 mL. / Max: 4500 Rpm. / RCF Max: 3024 xg / 30°.
1461-1743C		Angular 12 x 25/30 mL. / Max: 4500 Rpm. / RCF Max: 2716 xg / 30°.
1461-1746C		Angular 6 x 50 mL. / Max: 6000 Rpm. / RCF Max: 4427 xg / 30°.
1461-3275C		Angular 8 x 50 mL. / Max: 5000 Rpm. / RCF Max: 3494 xg / 30°. Biological security.
1461-1213		Angular 8 x 50 mL. / Max: 5000 Rpm. / RCF Max: 3494 xg / 30°.
1461-1744C		Angular 10 x 50 mL. / Max: 4500 Rpm. / RCF Max: 2830 xg / 30°.

Ref	Accessories	
1461-1718C		Angular 4 x 100 mL. / Max: 5000 Rpm. / RCF Max: 3158 xg / 30°.
1461-2200		Swing out 4 x 100 mL. / Max: 4000 Rpm. / RCF Max: 2469 xg.
1461-2200C		Swing out 4 x 100 mL. / Max: 4000 Rpm. / RCF Max: 2469 xg.
1461-2193C		Swing out 8 x 15 mL. / Max: 4000 Rpm. / RCF Max: 2504 xg.
1461-3215		Swing out 8 x 15 mL. Falcon / Max: 4000 Rpm. / RCF Max: 2469 xg.
1461-3113		Swing out 48 x 5 mL. / Max: 4000 Rpm. / RCF Max: 2469 xg.
1461-2218C		Swing out 2 microtiter plates / Max: 3000 Rpm. / RCF Max: 916 xg.
1461-2300		Horizontal 24 x 75mm. / Max: 13000 Rpm. / RCF Max: 16816 xg. Hematocrite.

Más centrifugas



Other centrifuges



BUNSEN SA. C/ Rio Alberche, 30 - 28970 Humanes de Madrid, Madrid, SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es

