



www.bunsen.es



ANALOG HEATING MANTLE

Ref. 1702

GENERAL CHARACTERISTICS

- * Electronic power regulator
- * Illuminated operating lights.
- * ABS case.
- * Excellent thermal insulation with fibroceraamic.

Manufactured according to EC directives.

Ref	Capacity	Dimensions Ø x H	Useful measures Ø x H	Max. temp	Power (W)
1702	500ml.	Ø230 x 155mm.	Ø115 x 70mm.	380°C	230



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fa. +34 916 128 254
info@bunsen.es - www.bunsen.es

