



www.bunsen.es



MINI-ORBITAL SHAKER

Ref. 1420

GENERAL CHARACTERISTICS

- * Shaker with vibratory movement with universal fixing system for different containers.
 - * Metal housing painted in epoxy.
 - * Electronic speed adjustment.
 - * Great stability at maximum shaking.
- * 60 minutes timer with a continuous switch function.
 - * Excentricity 6mm.
 - * Adjustable speed 100-500 Rpm.

Manufactured according to CE directives.




BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es



Ref	Vibration orbit	Consumption	Dimensions (W x H x D, mm)	Weight (Kg)	Max. load (Kg)	Useful platform (mm)
1420	6mm.	20 W.	230 x 115 x 295	6	2.5	220 x 220

Ref	Accessories
1450-1	25/50 mL. Flask adaptor
1450-2	100 mL. Flask adaptor
1450-3	250 mL. Flask adaptor
1450-4	500 mL. Flask adaptor
1450-5	1000 mL. Flask adaptor



BUNSEN SA. C/ Río Alberche, 30, - 28970 - Humanes de Madrid, Madrid. SPAIN
 Tlf. +34 916 113 584 - Fax. +34 916 128 254
info@bunsen.es - www.bunsen.es

